

WHY DIGITAL?

Digital sound processing offers many important benefits to improve your hearing experience. A key advantage of digital technology is that it allows your hearing healthcare professional to fit your hearing aid more precisely to your individual hearing loss and lifestyle.

WHY UNISON?

- ▶ Clean, crisp, comfortable sound
- ▶ Quick and accurate customization to your needs
- ▶ Helps maximize speech intelligibility in various listening situations

ABOUT UNITRON HEARING

Unison is part of the Unitron Hearing family of high-quality digital hearing aid solutions. For 40 years, Unitron Hearing has been committed to bettering the lives of the hearing impaired. Unitron Hearing is headquartered in Kitchener, Ontario, Canada, with operations in Canada, the United States and Germany, and serves customers in more than 50 countries around the world.

www.unitronhearing.com



unison



THE ADVANTAGES
OF DIGITAL
ARE NOW WITHIN
EVERYONE'S REACH





COMFORTABLE NATURAL SOUND

Unison offers many standard features and options that help ensure comfortable, natural sound, including:

- ▶ Quiet Mode Expansion – helps provide cleaner sound by reducing perceived noise within the hearing aid, caused by microphone noise or soft background sounds
- ▶ Feedback Manager – reduces the potential for uncomfortable feedback
- ▶ Telecoil – improves listening on the phone
- ▶ Manual Volume Control – provides full control over volume

The top section of the advertisement features a dark blue background. On the left, there are white, concentric ripples in water, suggesting sound waves. On the right, the word "unison" is written in a white, lowercase, sans-serif font, with a small white arc above the letter 'i'. Below the text, a row of five hearing aid shells is shown against the blue background. From left to right, they are: a light brown shell, a light brown shell, a light brown shell, a light brown shell, and a dark blue shell.

unison

JUST WHAT YOU WANT TO HEAR

A common complaint with hearing aids is the inability to understand speech in noisy environments. Unison's Vector Directional System addresses this by connecting you to the sounds in front of you and suppressing sounds from other directions. The Vector Directional System can be programmed into any of your listening programs and can even be implemented in smaller hearing aid styles.

LIVE IN UNISON

Unison is available in the full range of styles, sizes and colours to match your skin tone. Certain Unison features are available in specific shell styles only, so you should consult with your hearing healthcare professional to choose the one that best meets your needs.